



IntraOral CART

SOLUTIONS






Welcome to Medical Cart Company's Beaver range,
We created this range for the needs of the dental market.
You can find stands from the simplest use to real powered carts.



MAIN FEATURES




EASY TO CLEAN

Clean surfaces
without dirt traps.
One wipe cleaning.



ERGONOMIC

Specially designed
for dental procedures
and surgeries.




QUICK SETUP

Available assembled
and flat-packed,
easy to assemble.



SMALL FOOTPRINT

Developed for a small
area, fits in small
rooms too.



8 YEARS WARRANTY

Eight-year warranty
for the care-free
treatment.




LOCKING WHEELS

Simple locking
mechanism operated
by foot.

BEAVER 1

BV1


HEIGHT ADJUSTABLE
WORK STATION



BEAVER 1 PLUS

BV1-6-10


HEIGHT ADJUSTABLE
WORK STATION
WITH COMPUTER MOUNT



BEAVER 2

BV2


AIO COMPUTER
WORK STAND



BEAVER 2 PLUS

BV2-6-4-10


WORK STATION
WITH COMPUTER MOUNT



BEAVER 3

BV3-6-1


ALL-IN-ONE COMPUTER
WORK STATION



BEAVER 4

BV4-4-3


LAPTOP CART



BEAVER 5

2311-0618


SMALL DEVICE CART



BEAVER 6

BV6-4-9

PIEZOSURGERY CART




BEAVER 1

BV1
HEIGHT ADJUSTABLE WORKSTATION

-  *Ergonomic*
-  *Easy to clean*
-  *Small footprint*
-  *8 years warranty*
-  *Quick setup*
-  *Locking wheels*


Datasheet:
Click or Scan



FEATURES

HEIGHT ADJUSTMENT

Easy adjustment with gas spring. Pull the lever upward to release the lock and set it at the desired height.



INTRAORAL SCANNER HOLDER (OPTIONAL)

We installed a holder on the display mount for your IntraOral scanner.



DISPLAY MOUNT


Tiltable and rotatable VESA mount. For displays up to 24".

BEAVER 1 PLUS

BV1-6-10
HEIGHT ADJUSTABLE WORKSTATION
WITH COMPUTER MOUNT

-  *Ergonomic*
-  *Easy to clean*
-  *Small footprint*
-  *8 years warranty*
-  *Quick setup*
-  *Locking wheels*

Datasheet:
Click or Scan



FEATURES


ALL THE FEATURES OF
BEAVER1, PLUS:

KEYBOARD TRAY AND MOUSE PAD

Comfortable keyboard tray with a left or right hand mouse pad included. Enough space to use the mouse without limitations.

POWER SUPPLY HOLDER

You can place power supplies and chargers for your devices in the PSU holder.



COMPUTER MOUNT WITH POWER STRIP

Safe and secure mount for Tower/ Desktop PC. The straps keep the PC from falling.

You can power every device you need with the four-plug power strip.



BEAVER 2

BV2 ALL-IN-ON COMPUTER WORKSTAND



Ergonomic



Easy to clean



Small footprint



8 years warranty




Quick setup



Locking wheels

Datasheet:
Click or Scan




FEATURES

TILT/ROTATE DISPLAY MOUNT

With this solution you can adjust the screen freely. UP to 24" displays.

Tilt: +/-15°

Rotate: +/-90°



SCANNER HOLDER AND HANDLE

We created a solution for moving the stand easily. The scanner holder is on its front, just at the tip of your fingers.

KEYBOARD SHELF WITH PSU HOLDER

We created a big and comfortable shelf to place the keyboard with a mouse. In the PSU holder, you can place power supplies and chargers for your devices or even a 15" laptop.



BEAVER 2 PLUS

BV2-6-4-10

COMPUTER MOUNT WORKSTAND

-  *Ergonomic*
-  *Easy to clean*
-  *Small footprint*
-  *8 years warranty*
-  *Quick setup*
-  *Locking wheels*

Datasheet:
Click or Scan



FEATURES

ALL THE FEATURES OF
BEAVER2, PLUS:

COMPUTER MOUNT

Safe and secure mount for Tower/
Desktop PC.
The straps keep the PC from falling.



POWER STRIP

You can power every device you
need with the four-plug power
strip.



BEAVER 3

BV3-6-1
ALL-IN-ONE COMPUTER
WORKSTAND

 *Ergonomic*

 *Easy to clean*


 *Small footprint*

 *8 years warranty*

 *Quick setup*

 *Locking wheels*

Datasheet:
Click or Scan



FEATURES

TILT/ROTATE DISPLAY MOUNT

With this solution you can adjust the screen freely. UP to 24" displays.


Tilt: +/-15°

Rotate: +/-90°




SCANNER HOLDER AND HANDLE

We created a solution for moving the stand easily. The scanner holder is on its front, just at the tip of your fingers.



INTERNAL CABLE MANAGEMENT, POWER SUPPLY HOLDER

Hide all cables inside the column body. We added several holes for cables to enter and exit. Storage for the required power supplies.



BEAVER 4

BV4-4-3

LAPTOP CART

Suitable for 13" to 17" laptops.



Ergonomic



Easy to clean



Small footprint



8 years warranty



Quick setup



Locking wheels

Datasheet:
Click or Scan



FEATURES

LEFT-RIGHT HAND MOUSEPAD

Sliding mouse tray can be hidden under the worktop. Comfortable use for both left and right hand.



INTRAORAL SCANNER HOLDER AND SHELF

Small and clever solutions for intraoral scanner. We added a shelf for additional devices.

LAPTOP SECURITY

Optional laptop security to keep the laptop safe from theft and accidental falling.

The protector can be easily mounted on the worktop with knobs.



BEAVER 5

2311-0618

SMALL DEVICE CART

 *Ergonomic*

 *Easy to clean*

 *Small footprint*

 *8 years warranty*

 *Quick setup*

 *Locking wheels*

Datasheet:
Click or Scan



FEATURES

WORKTOP WITH HANDLE

Basic worktop for general use. You can place any instruments on it. We added a handle to make it easy to move.



SHELF

Height adjustable shelf to manage comfortable access. By loosening the screws one can adjust the height of the shelf. Folded up edges prevent objects to fall off.



BEAVER 6

BV6-4-9
PIEZOSURGERY CART

 *Ergonomic*

 *Easy to clean*

 *Small footprint*

 *8 years warranty*

 *Quick setup*

 *Locking wheels*



Datasheet:
Click or Scan



FEATURES

MULTIPLE SHELVES WITH POWER SUPPLY HOLDER

We created large shelves with storage under, where you can place your devices. There's enough space for piezosurgery tools without limitations.



HANDLE

We added a handle to the cart to let it move it easier.

CABLE MANAGEMENT

Multiple holes allow you to run the cables in the column, down to the bottom of the cart.

POWER STRIP

A four-plug powerstrip located on the bottom of the cart to power all you devices needed. Holes added for cable management.




MATERIAL AND PAINTING

We use aluminium as the main material for our carts. It's lightweight, strong and durable. After cutting and milling, the pieces get a powder coating for a better durability and ease of cleaning.

We paint the carts **RAL 9016 'Traffic white'** by default.

However we chose a pleasant color for the Beaver range:

TIGER Drylac® 68/50396



Custom colors are also available.
AVAILABLE FOR ALL CARTS



Medical Cart C O M P A N Y



www.dental-carts.com

Medical Cart Company Kft. - 2051 Hungary Biatorbágy, Budai út 16.
Phone: +36 23 920 790 - Email: info@medicalcartcompany.hu - www.medicalcartcompany.hu