

CHA6-0

Height adjustable cart
for general purposes

Chameleon 6 Range

Medical Cart
COMPANY



Ergonomic



Easy to clean



Multi purpose




5 years warranty



Quick setup




FEATURES



Built in channels to fix
accessories



Lockable wheels



Height adjustment



Designed and manufactured in
Hungary, products are typically
100% recyclable and upgradable
in the future. Raw materials are
sourced locally or in the EU.

Cart Properties


Dimensions and mechanical data




More info
on website




Accessories
available




Worksurface



Top



Front



Side

Mechanical data

Height	Min: 775 mm Max: 1140 mm
Wheelbase dimension	746 mm (diameter)
Work surface dimensions	Width: 620 mm Height: 560 mm
Wheels/locking	4 lockable, 100mm diameter
Loading capacity	35 kg

Construction

Main material	Aluminium, PMMA
---------------	-----------------

Color

RAL 9016	
Custom color	

Certificates

Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant