



MEDICAL  
CART COMPANY

Height adjustable cart  
with tablet mount

## Chameleon 6 Range



*Ergonomic*



*Easy to clean*



*Multi purpose*




*Lifetime warranty*



*Quick setup*




## FEATURES



*Height adjustment and  
channels to fix accessories*



*Lockable wheels*



*Tablet mount*



*Designed and manufactured in  
Hungary, products are typically  
100% recyclable and upgradable  
in the future. Raw materials are  
sourced locally or in the EU.*

# Cart Properties

## Dimensions and mechanical data






More info  
on website



Accessories  
available



Tablet Mount  
Datasheet



Mechanical data	
Work surface height	Min: 775 mm Max: 1135 mm
Wheelbase dimension	746 mm (diameter)
Work surface dimensions	Width: 620 mm Height: 560 mm
Wheels/locking	2 lockable, 100mm diameter
Loading capacity	35 kg
Construction	
Main material	Aluminium, PMMA
Color	
RAL 9016	
Custom color	
Certificates	
Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant