



MEDICAL
CART COMPANY

Height adjustable cart with
keyboard tray and display mount up to 24".

Chameleon 6 Range



Ergonomic



Easy to clean



Multi purpose




Lifetime warranty




Quick setup




FEATURES



*Height adjustment and
channels to fix accessories*



*Keyboard tray
with mouse pad*



*Tiltable, rotatable, height
adjustable Display Mount*



*Designed and manufactured in
Hungary, products are typically
100% recyclable and upgradable
in the future. Raw materials are
sourced locally or in the EU.*

Cart Properties

Dimensions and mechanical data



More info
on website






Accessories
available



Display Mount
Datasheet



Keyboard Tray
Datasheet



Mechanical data

Height	Min: 775 mm Max: 1135 mm
Wheelbase dimension	746 mm (diameter)
Work surface dimensions	Width: 620 mm Height: 560 mm
Wheels/locking	2 lockable, 100mm diameter
Loading capacity	35 kg
VESA Mount hole pattern	75x75 mm; 100x100 mm
VESA Mount load capacity	6-12 kg up to 24"

Construction

Main material	Aluminium, PMMA
---------------	-----------------

Color

RAL 9016	
Custom color	

Certificates

Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant