



RHI-0

Height adjustable cart
for heavy duty purposes

Rhino Range



Ergonomic



Easy to clean



Multi purpose



Lifetime warranty




Quick setup




FEATURES



Designed and manufactured in Hungary, products are typically 100% recyclable and upgradable in the future. Raw materials are sourced locally or in the EU.



Lockable wheels



*Height adjustment
and column channels
to fix accessories*

Cart Properties



Dimensions and mechanical data



More info
on website



Accessories
Available



Mechanical data

Height	Min.: 865 mm Max.: 1200 mm
Wheelbase dimension	750 mm (diameter)
Work surface dimensions	Width: 620 mm Height: 560 mm
Wheels/locking	2 lockable, 125mm diameter
Loading capacity	80 kg

Construction

Main material	Aluminium
---------------	-----------

Color

RAL 9016	
Custom color	

Certificates

Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant