

# CM2 - Camera Mount with Display Mount

With adjustable VES2 VESA head.

MEDICAL  
CART COMPANY

Product range:  



## **Ergonomic**

*Tiltable, rotatable,  
height adjustable*



## **Easy to clean**

*Clean surfaces with  
minimal dirt traps*



## **Durable**

*Lifetime warranty*



## **Quick setup**

*Supplied fully assembled  
with minimal packaging*



## **Attachable display**

*Up to 24 inches*




## **Load capacity**


*From 6 kg to 12 kg*




## FEATURES



*Tiltable*



*Rotatable*



*Height adjustable camera  
and display mount*

# Display Mount Properties

## Dimensions and mechanical data




More info  
on website


### Technical attributes:

- Fully recyclable
- Lifetime Warranty
- Made in Hungary


### Product range:



Front



Top



Front

### Cable management:



Tangle-free solution that helps to run the cables in the column.

Mechanical data	
Camera mount height	550-850 mm
Column height	500 mm
Hole pattern	75x75mm; 100x100 mm
Net. Weight	2.5 kg
Load capacity	6-12 kg
Attachable display	Up to 24"
Construction	
Material	Steel, Aluminium
Color	
RAL 9016	
Custom color	
Certificates	
Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant