

Computer holder 1

CPU1

Product range: 


MEDICAL
CART COMPANY

Holder Properties

Dimensions and mechanical data



More info
on website



Durable
Lifetime warranty



Quick setup
Can be added to the cart quick
and easily



Easy to clean
Clean surfaces with
minimal dirt traps




Ergonomic
Easy to install any size desktop
computer on it

FEATURES




Can be mounted on any column




Can be adjusted for bigger sizes

Technical attributes:

- Fully recyclable
- Lifetime Warranty
- Made in Hungary



Side



Bottom

| Mechanical data | |
|--------------------------|---|
| Width | 250 mm |
| Height | 221 mm |
| Depth | 47 mm - 97 mm |
| Max. weight load | 4 kg |
| Construction | |
| Main material | Steel |
| Color | |
| RAL 9016 |  |
| Custom color |  |
| Certificates | |
| Certificates conforms to | IEC60601-1:2012 (Edition 3) |
| CE | Rohs compliant |