


Computer holder 3

CPU3

MEDICAL
CART COMPANY

Product range:    



Durable

Lifetime warranty



Quick setup

Can be added to the cart quick and easily



Easy to clean

Clean surfaces with minimal dirt traps



Ergonomic

Easy to install any size desktop computer on it

FEATURES



Can be mounted on any column



Adjusted for different sizes

Holder Properties

Dimensions and mechanical data




More info
on website


Technical attributes:


- Fully recyclable
- Lifetime Warranty
- Made in Hungary

Side



Bottom



Mechanical data	
Width	250 mm
Height	221 mm
Depth	130 mm - 214 mm
Max. weight load	5 kg
Construction	
Main material	Steel
Color	
RAL 9016	
Custom color	
Certificates	
Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant