

# DM5 - Display Mount

With VES1 VESA head.

MEDICAL  
CART COMPANY

Product range:  



**Ergonomic**



**Easy to clean**

Clean surfaces with minimal dirt traps



**Durable**

Lifetime warranty



**Quick setup**

Supplied fully assembled with minimal packaging



**Attachable display**

Up to 27 inches




**Load capacity**


Up to 15 kg



## FEATURES



**Teardrop holes to hang display quickly**



**Height adjustable**

# Display Mount Properties

## Dimensions and mechanical data






More info  
on website

### Technical attributes:

- Fully recyclable
- Lifetime Warranty
- Made in Hungary

### Product range:



Mechanical data	
Total height	430 mm
Column height	Standard: 400 mm; Optional: 500 mm
Hole pattern	75x75 mm; 100x100 mm
Net. Weight	1 kg
Load capacity	Up to 15 kg
Attachable display	Up to 27"
Construction	
Material	Steel, Aluminium
Color	
RAL 9016	
Custom color	
Certificates	
Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant