

DM6 - Display Mount

With VES1 VESA head.

MEDICAL
CART COMPANY

Product range:  



Ergonomic



Easy to clean

Clean surfaces with minimal dirt traps



Durable

Lifetime warranty



Quick setup

Supplied fully assembled with minimal packaging



Attachable display

Up to 27 inches



Load capacity


Up to 20 kg



FEATURES



Teardrop holes to hang display quickly



Height adjustable on both sides

Display Mount Properties

Dimensions and mechanical data




More info
on website


Technical attributes:

- Fully recyclable
- Lifetime Warranty
- Made in Hungary


Product range:



Front



Top



Front

Mechanical data

Total height	430 mm
Column height	Standard: 400 mm; Optional: 500 mm
Hole pattern	75x75 mm; 100x100 mm
Net. Weight	1.2 kg
Load capacity	Max. 20 kg
Attachable display	Max. 27" on each VESA head.

Construction

Material	Steel, Aluminium
----------	------------------

Color

RAL 9016	
Custom color	

Certificates

Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant