

BUDGIE CART SOLUTIONS 2025

Medical Cart
COMPANY



GENERAL FEATURES



EASY TO CLEAN SURFACES

Clean surfaces without dirt traps.
One wipe cleaning.



QUICK SETUP

Available assembled or
flat-packed, easy to assemble.



8 YEARS WARRANTY

Eight-year warranty for
the care-free treatment.



ERGONOMIC

Designed for easy
and comfortable use.



SMALL FOOTPRINT

Developed to fit narrow doors,
ease of movement.



LOCKABLE WHEELS

Simple locking mechanism
operated by foot.



CORE PRODUCT

BUG-0

BASIC BUDGIE CART

- 100mm locking wheels
- B1 base - 25kg max loading capacity
- C2 lightweight column

Datasheet
Scan the QR



PRODUCT WITH ACCESSORIES

BUG-6

BUDGIE CART with Quick release mount

- Quick release mount
- Handle
- DIN rail
- Drip pole
- Humidifier mount
- Basket
- PSU Holder



BUG-0

Basic lightweight cart.



Datasheet



Scan the QR

BUG-6

Lightweight cart with Quick release mount.



Datasheet



Scan the QR

BUG-7

Lightweight cart with Drip pole.



Datasheet



Scan the QR

WORKTOP OPTIONS

WS2

Material:
Aluminium


Dimensions:
W: 400mm
D: 450mm
THK: 4mm



WS3

Material:
Acryl

Dimensions:
W: 430mm
D: 450mm
THK: 8mm



WS12

Material:
Acryl

Dimensions:
W: 560mm
D: 485mm
THK: 8mm



WS13

Material:
Acryl

Dimensions:
W: 460mm
D: 340mm
THK: 8mm



WS-US3

Material:
Acryl


Dimensions:
W: 350mm
D: 350mm
THK: 8mm




MNT-1

Material:
Aluminium

Dimensions:
W: 430mm
D: 565mm
THK: 6mm



**Custom
design
on
request**



ACCESSORIES

DISPLAY MOUNT




VES1

Datasheet
Scan the QR



TABLET MOUNT




TM2

Datasheet
Scan the QR



QUICK RELEASE



QRA 1/2/3/4/5/6

Datasheet
Scan the QR



BASKET



ST 1/6/8

Datasheet
Scan the QR



GAS BOTTLE HOLDER




GBH3

Datasheet
Scan the QR



SHELVES




S1, S2, S3, S11

Datasheet
Scan the QR



SCANNER HOLDER



SH 1/2/3/4

Datasheet
Scan the QR



PSU HOLDER



PSH 1/2/4/5/7/9

Datasheet
Scan the QR



POWER STRIP



PDU 3/4/5/8

Datasheet
Scan the QR



COMPUTER MOUNT



CPU 1/2

Datasheet
Scan the QR



DIN RAIL



DR2, DR3

Datasheet
Scan the QR



DRIP POLE



IV1-DRC

Datasheet
Scan the QR



RUBBER HOLDER




SR1, SR2, SR3

Datasheet

Scan the QR



HUMIDIFIER MOUNT




AM1

Datasheet

Scan the QR



POAG STUD



POW-A2-15/25

Datasheet

Scan the QR



CABLE HOOK



CH 1/2/3/4

Datasheet

Scan the QR



MOUSE HOLDER



MH1


Datasheet

Scan the QR



CUSTOM COLORS

We offer the choice of custom colors for the aluminium components from **Tiger 3D metallic range**.



Visit the Tiger webshop and select the color you desire:

<https://www.tiger-coatings.com>

Medical Cart C O M P A N Y



Medical Cart Company Kft.

2051 Hungary

Biatorbágy, Budai út 16.

Phone: +36 23 920 790

Email: info@medicalcartcompany.hu

www.medicalcartcompany.hu