

CHAMELEON CART SOLUTIONS 2025

Medical Cart
COMPANY



GENERAL FEATURES



EASY TO CLEAN SURFACES

Clean surfaces without dirt traps.
One wipe cleaning.



QUICK SETUP

Available assembled or
flat-packed, easy to assemble.



8 YEARS WARRANTY

Eight-year warranty for
the care-free treatment.



ERGONOMIC

Designed for easy
and comfortable use.



SMALL FOOTPRINT

Developed to fit narrow doors,
ease of movement.



LOCKABLE WHEELS

Simple locking mechanism
operated by foot.



Datasheet
Scan the QR



CORE PRODUCTS

CHA4-0

BASIC CHAMELEON 4 CART

- B1 base - 25kg max loading capacity, 100mm wheels
- WS6 work surface




CHA4-1
DISPLAY MOUNT



CHA4-2
TABLET MOUNT



CHA4-3
LAPTOP MOUNT



CHA4-4
DISPLAY MOUNT
+
KEYBOARD TRAY


CHA6-0

BASIC CHAMELEON 6 CART

- B3 base - 35kg max loading capacity, 100mm wheels
- WS9 work surface




CHA6-1
DISPLAY MOUNT



CHA6-2
TABLET MOUNT



CHA6-3
LAPTOP MOUNT



CHA6-4
DISPLAY MOUNT
+
KEYBOARD TRAY

CHA4-3-B

Budget laptop cart.



Datasheet



Scan the QR

CHA4-8

Large work surface cart.



Datasheet



Scan the QR

CHA6-4-CPU

Complex workstation with CPU mount, Display mount and Keyboard tray.



Datasheet




Scan the QR

WORKTOP OPTIONS

WS3

Material:
Acryl

Dimensions:
W: 430mm
D: 450mm
THK: 8mm



WS6

Material:
Acryl


Dimensions:
W: 540mm
D: 485mm
THK: 8mm



WS7

Material:
Acryl


Dimensions:
W: 620mm
D: 560mm
THK: 8mm



WS9

Material:
Acryl

Dimensions:
W: 620mm
D: 560mm
THK: 8mm



WS13

Material:
Acryl


Dimensions:
W: 460mm
D: 340mm
THK: 8mm



WS15

Material:
Acryl

Dimensions:
W: 700mm
D: 500mm
THK: 8mm




1-0770

Material:
Aluminium

Dimensions:
W: 500mm
D: 560mm
THK: 8mm



**Custom
design
on
request**



ACCESSORIES

DISPLAY MOUNT




DM 2/5/6/7/8/9

Datasheet
Scan the QR



TABLET MOUNT



TM2

Datasheet
Scan the QR



CAMERA MOUNT



CM1, CM2

Datasheet
Scan the QR



STORAGE




ST 1/2/3/5/6/7/8

Datasheet
Scan the QR



SHELVES




S 1/2/3/4/6/8/10/11

Datasheet
Scan the QR



DRAWERS




D 1/2/3/4/5/6/7/8/9/10/11

Datasheet
Scan the QR



SCANNER HOLDER




SH 1/2/3/4

Datasheet
Scan the QR



PSU HOLDER



PSH 1/2/4/5/7/8/9

Datasheet
Scan the QR



POWER STRIP



PDU 3/4/5/8

Datasheet
Scan the QR



COMPUTER MOUNT




CPU 1/2/3/4

Datasheet
Scan the QR



DIN RAIL




DR1, DR3

Datasheet
Scan the QR



POLE CLAMP




DRC 1/2

Datasheet
Scan the QR



ACCESSORIES

RUBBER HOLDER




SR1, SR2, SR3

Datasheet

Scan the QR



HUMIDIFIER MOUNT




AM1

Datasheet

Scan the QR



POAG STUD



POW-A2-15/25

Datasheet

Scan the QR



POAG STUDS



POW-A1

Datasheet

Scan the QR



CABLE HOOK



CH 1/2/3/4

Datasheet

Scan the QR



MOUSE HOLDER



MH1

Datasheet

Scan the QR



COUNTERWEIGHT



CBW-3, CBW-6


Datasheet

Scan the QR



CUSTOM COLORS

We offer the choice of custom colors for the aluminium components from **Tiger 3D metallic range**.



Visit the Tiger webshop and select the color you desire:

<https://www.tiger-coatings.com>

Medical Cart C O M P A N Y



Medical Cart Company Kft.

2051 Hungary

Biatorbágy, Budai út 16.

Phone: +36 23 920 790

Email: info@medicalcartcompany.hu

www.medicalcartcompany.hu