

RHINO CART SOLUTIONS 2025

 Medical Cart
COMPANY



GENERAL FEATURES



EASY TO CLEAN SURFACES
Clean surfaces without dirt traps.
One wipe cleaning.



QUICK SETUP
Available assembled or
flat-packed, easy to assemble.



8 YEARS WARRANTY
Eight-year warranty for
the care-free treatment.



ERGONOMIC
Designed for easy
and comfortable use.



SMALL FOOTPRINT
Developed to fit narrow doors,
ease of movement.



LOCKABLE WHEELS
Simple locking mechanism
operated by foot.

CORE PRODUCTS

RHI-0

BASIC RHINO CART

- B6 base - 80kg max loading capacity, 125mm wheels
- WS9 work surface



Datasheet
Scan the QR



RHI-1 DISPLAY MOUNT



RHI-4 DISPLAY MOUNT + KEYBOARD TRAY



RHI-4-POW

Powered cart with battery pack, large PSU holder and power strip.



Datasheet



Scan the QR


DIGITAL FEEDBACK PANEL

Placed on the work surface.
Service related conditions shown on the display.




POWERED BASE

Battery pack can be mounted from below.



ISOLATION TRANSFORMER HOLDER

Simple holder option for the most common Isolation transformers on the market.



DRAWERS

Multiple drawer options:
Small and large, DIN rail mounting,
key lock, push-to-open.

Datasheet



Scan the QR




CABLE MANAGEMENT

Folded sheet component can be
mounted on the back of the column
to keep the cables collected.



LARGE LOWER SHELF

For CPU, printer or larger devices.



WORKTOP OPTIONS

1-0770

Rear handle

Material:
Aluminium

Dimensions:

W: 500mm
D: 560mm
THK: 8mm



1-1415

Front handle


Material:
Aluminium

Dimensions:

W: 500mm
D: 560mm
THK: 8mm



**Custom
design
on
request**



ACCESSORIES

DISPLAY MOUNT




DM 2/5/6/7/8/9

Datasheet
Scan the QR



TABLET MOUNT



TM2

Datasheet
Scan the QR



CAMERA MOUNT



CM1, CM2

Datasheet
Scan the QR



STORAGE




ST 1/2/3/5/6/7/8

Datasheet
Scan the QR



SHELVES



S 1/2/3/4/6/8/10/11

Datasheet
Scan the QR



DRAWERS




D 1/2/3/4/5/6/7/8/9/10/11

Datasheet
Scan the QR



SCANNER HOLDER




SH 1/2/3/4

Datasheet
Scan the QR



PSU HOLDER



PSH 1/2/4/5/7/8/9

Datasheet
Scan the QR



POWER STRIP



PDU 3/4/5/8

Datasheet
Scan the QR



COMPUTER MOUNT




CPU 1/2/3/4

Datasheet
Scan the QR



DIN RAIL




DR1, DR3

Datasheet
Scan the QR



POLE CLAMP




DRC 1/2

Datasheet
Scan the QR



ACCESSORIES

RUBBER HOLDER



SR1, SR2, SR3

Datasheet

Scan the QR



HUMIDIFIER MOUNT




AM1

Datasheet

Scan the QR



POAG STUD



POW-A2-15/25

Datasheet

Scan the QR



CABLE HOOK



CH 1/2/3/4

Datasheet

Scan the QR



MOUSE HOLDER



MH1

Datasheet

Scan the QR



POAG STUDS



POW-A1


Datasheet

Scan the QR



CUSTOM COLORS

We offer the choice of custom colors for the aluminium components from **Tiger 3D metallic range**.



Visit the Tiger webshop and select the color you desire:

<https://www.tiger-coatings.com>

Medical Cart C O M P A N Y



Medical Cart Company Kft.

2051 Hungary

Biatorbágy, Budai út 16.

Phone: +36 23 920 790

Email: info@medicalcartcompany.hu

www.medicalcartcompany.hu