

MEDICAL CART COMPANY



The **Medical Cart Company Ltd** was born in the year 2016 and being a proud Hungarian company, all products are designed in house and material sourced and manufactured locally or within Europe. We have a total commitment to reducing our carbon footprint.



We use materials which are 100% recyclable and ensure that all products have the ability to be refurbished after many years of use and adapted to the changing needs of users with minimal effort.

How to choose your cart ?

We have designed a range of bases to suit the requirements of your equipment and the accessories required to support good use and function. All products are available with and without height adjustment .

So how to choose

- Weight of the equipment to accessories to be supported
- Working height required
- Size of work surface
- Location of accessories

We have a standard range of carts which will suit most needs or alternatively we can design the cart to your needs.

We understand, that a standard mounting solution may not your need, so we offer a fast and affordable customisation service as well as an exclusive for you and OEM design service.

We are here to be part of your team, not just your supplier !




Contents

- Budgie range 3
- Giraffe range 5
- Chameleon 4, 5, 6 range 7
- Rhino range 9
- Drawers 11
- Snake range 13
 - Wall Mount Solutions* 15
 - Circuit Arm* 16
- Sit/Stand wall Mount 17
- Desk Stands 18
- Components for Sale 19
- Accessories 21
- Power Solutions 23
- Power Pack 24
- Jungle Carts 25
- Custom Carts 26



Budgie range



Features

Elegant, strong, light weight and safe roll stand to suit the differing needs of your requirements. A choice of bases ensures that most applications and weight requirements can be served.



Ergonomic



Easy to Clean



Durable



Multi Purpose



Quick Setup



Color Configuration




B1

Max Load 36kg



B2

Max Load 22kg



Scan or Click QR code for datasheet

Popular Configurations



IV Pole Cart

Accessories: Basket, Drip Pole, Cable Hook, Handle, B2 Base



Tablet Mount Cart

Accessories: Tablet mount, Din Rail, Gel Holder, B2 Base



Quick Release Cart

Accessories: Quick Release Head, Din Rail, Handle and Storage bin, Storage Bin, B2 Base



Display Mount Cart

Accessories: Vesa Mount, Din Rail, B2 Base



Giraffe range



Features

Strong, light weight roll stand to suit the differing needs of your requirements. A choice of bases ensures that most applications and weight requirements can be served.



Ergonomic



Easy to Clean



Durable



Multi Purpose



Quick Setup



Color Configuration



B1

Max Load 36kg



B3

Max Load 65kg



B4

Max Load 85kg



B5


Max Load 40kg



B6

Max Load 80kg

Optional Vesa 75/100 Display Mount/ Tablet Mount



Optional ESD Wheels



Scan or Click QR code for datasheet

Popular Configurations



Quick Release Head Cart

Accessories: Quick Release Head, Din Rail, Drip Pole, Basket, B1 Base



Height Adjustable Stand upto 50" Displays

Accessories: Vesa Mount, Handle, Basket, B1 Base



Ventilator Cart

Accessories: Top Plate, Handle, Cabel Hook, Basket, Din Rail, B3 Base



Tablet Mount Cart

Accessories: Tablet mount, Din Rail, Gel holder, Handle, Storage Bin ,PSU holder, B1 Base



Display Mount Cart

Accessories: Display mount, Plastic Storage Bin, Shelf, Basket, PSU holder, Hook, B3 Base



Chameleon range



Info

- Chameleon 4** (B1 base which max load 36kg, size 616mm and small worksurface)
- Chameleon 5** (B5 base which max load 40kg, size 697mm and large worksurface)
- Chameleon 6** (B3 base which max load 65kg, size 746mm and large worksurface)


Ergonomic	Easy to Clean	Durable
Multi Purpose	Quick Setup	Color Configuration

Popular Configurations



Ultrasound Cart

Accessories: Storage Bin, Shelf, B3 Base, Aluminium Work Surface WS12 work surface



Camera Mount Cart

Accessories: Camera Mount, Keyboard tray, Din Rail Kit, Mouse Holder, Drawer, PSU Holder, B5 Base, WS11 Work Surface



Tablet Mount Cart

Accessories: Tablet Mount, PSU Holder, Drawers, B3 Base, WS11 Work Surface




B1
Max Load 36kg



B3
Max Load 65kg



B5
Max Load 40kg



Laptop Mount Cart

Accessories: Laptop Mount, Keyboard Tray, Din Rail Kit, Mouse Holder, 2nd Work Surface, PSU Holder, B1 Base, WS5 Work Surface



ECG Cart

Accessories: Module Mount, Din Rail Kit, Gel Holder, B1 Base, WS5 Work Surface



Display Mount Cart

Accessories: Display Mount, Keyboard tray, Din Rail Kit, Gel Holder, Mouse Holder, Basket, PSU Holder, B1 Base, WS5 Work Surface



Scan or press QR code datasheet



Rhino range



Features

Strong to suit the differing needs of your requirements. A choice of bases ensures that most applications and weight requirements can be served.

- Ergonomic
- Autonomy
- Durable
- Multi Purpose
- Easy to Clean
- Color Configuration
- Quick Setup

Choice of Work Surface

Large Range of Accessories

Cable Management

Optional ESD Wheels

Height Adjustment

Range of Mounting Platforms



B6
Max Load 80kg




Scan or Click QR code for datasheet

Popular Configurations



Camera Mount Cart

Accessories: Camera Mount, Keyboard tray, Din Rail Kit, Mouse Holder, PSU Holder, B6 Base, WS11 Work Surface



Display Mount Cart

Accessories: Display Mount, Din Rail Kit, Mouse Holder, Gel Holder, PSU Holder, Shelf, B6 Base, WS11 Work Surface



Neuro Cart

Accessories: Din Rail Kit, Keyboard tray, Din Rail, Mouse Holder, PSU Holder, B6 Base, Aluminium Work Surface



Equipment Cart

Accessories: Display Mount, Keyboard tray, Din Rail, Mouse Holder, B6 Base, Aluminium Work Surface

Drawers




Info

The Drawers easy to clean and easy to assemble.
Only suitable for Chameleon 5,6 and Rhino

- 
Multi Purpose
- 
Quick Setup
- 
Easy to Clean
- 
ISO Tray Option


Small Drawer (80mm)



Large Drawer (160mm)




Small + Large




Scan or Click QR code for datasheet

ISO Tray 40x30cm



Tray for Partitons
300x400x50mm
*See related partitions fo tray



Tray
300x400x100mm



Optional Inserts

INDIVIDUAL MODULAR INSERTS




Individual modular inserts arranged
in a tray


2162205 + 2152206



2162205 Partition 70mm
for 2162206
2162206 Inserts, 170x75x50mm



2862751 + 2862753




Suspension rails for inserts 2162206
2862751 550mm, (fits 7 inserts)
2862753 650mm, (fits 8 insterts)




2162204 + 2162202




Inserts:
2162204 100x50x45mm
2162202 50x50x45mm



2162210



Inserts:
120x110x50mm



2162220 + 2162212



Inserts:
2162220 240x55x50mm
2162212 240x110x50mm




Partition option


PARTITION "SNAP" FEATURE



All partitions have a "snap"
feature between the ribs, for
defined division of the modules.



PARTITIONS FOR BASKETS, ABS





Snake range



Wall Mount For Display and Keyboard


Following the lead of the industry leaders we now have compatible solutions to fit into existing wall channels. The arms have 180 degrees rotation, hidden cable management. The arms are supplied assembled and very easy to use and clean.



WA-F-KD-E2
Display and Keyboard Mount
Double Extension



WA-H-KD-E1
Display and Keyboard Mount Height
adjustable, single Extension



WA-H-KD
Display and Keyboard Mount
Height Adjustable

Wall mount solutions


Single or double arms suitable for displays and medical equipment, hidden cable management and customisation is available.



WA-D



WA-F-D-E1
Extension



WA-F-D-E2
Extension



WA-K




WA-F-K-E1
Keyboard Single
Extension




WA-F-K-E2
Keyboard Double
Extension

Cable Management


Organised Cable management enables easier cleaning techniques.



Circuit Arm



A well designed heavy duty circuit arm ,designed with a constant friction hinge which can be adjusted to the correct force correct of a specific application. Custom lengths and fitting can be made and supplied.



Tube Holder

Allows for easy installation without damage to the tubes or cables to be supported.



Articulated Joint

Constant force joint with the ability to adjust the force enables a safe setup.




Mounting options

Standard Din Rail



Tube Holder




Sit/Stand wall Mount



Foldaway keyboard tray, allows busy areas to be used as a work station without consuming significant floor space.



Desk Stands


-  Height Adjustable Display Mount
-  Integrated Cable Management
-  Tilt and Rotation



Components Sales


Work surface

Our standard range of work surfaces can be fitted to any of our carts or alternatively fitted to your design. Through customers requests, these are the most requested sizes and functions of the work surface integrating handles, scanner holders and cable holders.




Aluminum Extrusions

6mtr lengths of aluminium profile, can be supplied as processed or raw.



Base

Bases Cast aluminum and can be supplied finished or raw and are available to purchase separately can be finished in any colour. Weight is defined with no height adjustment with loaded weight upto 1mtr from the floor.





Wheels


We offer 2 versions of Twin wheels , economy range and the premium range from Steincop plus ESD options. All, wheels are supplied with a M12 thread.



Accessories

	Display Mount DM1  
	Tablet Mount TM1  
	Laptop Mount LM1  
	Camera Mount CM1  
	Quick Release Head QR1  
	Drip Pole IV1  

	2nd Work Surface 40 2WS4  
	2nd Work Surface 50 2WS5  
	Keyboard Tray KB1  
	Keyboard Tray KB2  

	
	Storages  
	

	
	Shelves  
	
	

	CPU Holder CPH1  
	Mouse Holder MH1  
	Hook Cable Holder CH1  
	Scanner holders SH1, SH2  
	Gel Holders GH1, GH2  

	
	PSU Holders  
	
	Handle H1  
	Fisher and Paykel AM1  


	Side Din Rail Kit DR1  
	Top Plate with Din Rail 1911-9201  
	Din Rails DR2, DR3, DR4, DR5, DR6  
	Din Rail Mounts DRC1, DRC2  



Scan or Click QR code for datasheet

Power Solutions

Medical Isolation Transformers for added safety when using medical technology



REOMED 300VA A-50W Series Isolation Transformer



REOMED 600VA A-50W Series Isolation Transformer




REOMED 1000VA A-50W Series Isolation Transformer



Medical Power Strip



Potential Equalization Cable 3 m / 5 m



Potential Equalization Socket



Potential Equalization Pin 15 mm



Potential Equalization Pin 25 mm

PWP Power System



Power Pack information

Heavy Duty LiPo4 Power Pack suitable to fit any cart with built in charger and AC inverter allow easy installation of suitable equipment. Available in 512Wh or 240Wh sizes.



Autonomy



Easy to Clean




Durable



Multi Purpose



Ease of Use



Intelligent fuel gauge report run time and battery health. Datalogging and fleet management software is built into the product as standard.

LED Indicator

Easy to see, easy to understand battery charge indicator.




Service Related Conditions

Service related alarms, such as overload, over temperature, replace battery etc., will cause one of the RUI Service LEDs to flash.



Jungle Carts



We believe that patients and users should have the fun of alternative colours and artwork on the carts, we offer a range of standard designs or we can print and design to your requirement.



Custom Carts

3D CONCEPT

We of the concept and functionality provide 3D realization of your ideas to help you get a view of the concept. The functionality, usability and modularity of our products can be tested via the prototype designed in accordance with your needs.



PRODUCTION

After the testing phase based on the standards and quality expectations, production starts. Taking account of your budget constraints and constant target, our robust production capacities ensure best-in-class products and reliable delivery.



MEDICAL CART COMPANY

Medical Cart Company Ltd.
Hungary, 2051 Biatorbágy
Vendel park, Tormásrét utca 12.
Tel: +36 23 920 790
Tel: +36 30 841 4353

www.medicalcartcompany.hu

11 Okt. 2021 v 2.0

