

# **BUG-0-3-0038**

**MEDICAL**  
CART COMPANY

## **Accutronics VyairFabian Evo Cart** DIN Rails, Handle, Hook and Humidifier mount



*Ergonomic*



*Easy to clean*



*Multi purpose*



*Lifetime warranty*




*Quick setup*




## **FEATURES**



*Designed and manufactured in Hungary, products are typically 100% recyclable and upgradable in the future. Raw materials are sourced locally or in the EU.*



*Safe device mounting*



*Channels to fix accessories*



*Lockable wheels*

# Cart Properties

## Dimensions and mechanical data




More info  
on website




Accessories  
available




Handle  
Datasheet




Mounting plate




DIN Rail




Humidifier mount



Base



Front



Side

### Mechanical data

Top plate dimensions	232x220 mm
Height	955 mm
Wheelbase dimension	616 mm (diameter)
Wheel construction	Polyamide with Polyurethane coating
Wheels/locking	2 lockable, 100 mm diameter
Loading capacity	25 kg

### Construction

Main material	Aluminium
---------------	-----------

### Color

RAL 9016	
Custom color	

### Certificates

Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant