

Ultrasound & EKG Cart

MEDICAL
CART COMPANY

Budgie Range

BUG-0-WS12



Ergonomic



Easy to clean



Multi purpose



Lifetime warranty




Quick setup




FEATURES



Designed and manufactured in Hungary, products are typically 100% recyclable and upgradable in the future. Raw materials are sourced locally or in the EU.



*Channels to
fix accessories*



*Optional rubber inserts
for accessories*




Lockable wheels

Cart Properties


Dimensions and mechanical data




More info
on website




Work surface



Base



Front



Side

Mechanical data

Height	1030 mm
Wheelbase dimension	616 mm (diameter)
Wheel construction	Polyamide with Polyurethane tread
Wheels/locking	2 lockable, 100 mm diameter
Loading capacity	36 kg
Work surface size	560x485 mm

Construction

Main material	Aluminium
---------------	-----------

Color

RAL 9016	
Custom color	

Certificates

Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant