



VESA Mount Budgie Cart



- Ergonomic
- Easy to clean
- Multi purpose
- Lifetime warranty
- Quick setup




FEATURES



Designed and manufactured in Hungary, products are typically 100% recyclable and upgradable in the future. Raw materials are sourced locally or in the EU.



Channels to fix accessories



VESA mount



Lockable wheels

Cart Properties

Dimensions and mechanical data




More info on website




Accessories available




VESA Mount Datasheet




VESA Head up to 27"



Base



Front



Side

Mechanical data

Height	1500 mm
Wheelbase dimension	616 mm (diameter)
Wheel construction	Polyamide with Polyurethane tread
Wheels/locking	2 lockable, 100 mm diameter
Loading capacity	36 kg
VESA Mount hole pattern	75x75 mm, 100x100 mm
Display size	Up to 27", max. 15 kg

Construction

Main material	Aluminium
---------------	-----------

Color

RAL 9016	
Custom color	

Certificates

Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant