

Drip Pole Cart

MEDICAL
CART COMPANY

Budgie Range

BUG-7



Ergonomic



Easy to clean




Multi purpose



Lifetime warranty




Quick setup



FEATURES



Designed and manufactured in Hungary, products are typically 100% recyclable and upgradable in the future. Raw materials are sourced locally or in the EU.



*Channels to
fix accessories*



Drip pole



Lockable wheels


Cart Properties

Dimensions and mechanical data




More info
on website


Load capacity of hooks:
1 kg!



Base



Front



Side

Mechanical data

Height	Min.: 1350 mm; Max.: 2200 mm
Wheelbase dimension	616 mm (diameter)
Wheel construction	Polyamide with Polyurethane tread
Wheels/locking	2 lockable, 100 mm diameter
Loading capacity	36 kg
Hook maximum load	1 kg/hook

Construction

Main material	Aluminium
---------------	-----------

Color

RAL 9016	
Custom color	

Certificates

Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant