




## Drip pole Budgie Cart



- Ergonomic**
- Easy to clean**
- Multi purpose**
- Lifetime warranty**
- Quick setup**

### FEATURES



Designed and manufactured in Hungary, products are typically 100% recyclable and upgradable in the future. Raw materials are sourced locally or in the EU.



Channels to fix accessories



Drip pole



Lockable wheels

# Cart Properties

## Dimensions and mechanical data




More info on website


Accessories available

VESA Mount Datasheet


**Load capacity of hooks:  
1 kg!**



Base



Front



Side

Mechanical data	
Height	Min.: 1350 mm; Max.: 2200 mm
Wheelbase dimension	616 mm (diameter)
Wheel construction	Polyamide with Polyurethane tread
Wheels/locking	2 lockable, 100 mm diameter
Loading capacity	36 kg
Hook maximum load	1 kg/hook
Construction	
Main material	Aluminium
Color	
RAL 9016	
Custom color	
Certificates	
Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant