



MEDICAL  
CART COMPANY

Height adjustable cart  
with laptop mount

## Chameleon 4 Range



*Ergonomic*



*Easy to clean*



*Multi purpose*




*Lifetime warranty*



*Quick setup*




## FEATURES



*Height adjustment and  
channels to fix accessories*



*Lockable wheels*



*Laptop protection up to 17"*




*Designed and manufactured in  
Hungary, products are typically  
100% recyclable and upgradable  
in the future. Raw materials are  
sourced locally or in the EU.*

# Cart Properties


## Dimensions and mechanical data



More info  
on website



Work surface




Laptop stand




Accessories  
available




Laptop Mount  
Datasheet



Top



Front



Side

### Mechanical data

|                         |                                          |
|-------------------------|------------------------------------------|
| Work surface height     | Min: 775 mm Max: 1135 mm                 |
| Wheelbase dimension     | 616 mm (diameter)                        |
| Work surface dimensions | Width: 540 mm Height: 485 mm             |
| Wheels/locking          | 2 lockable, 100 mm diameter              |
| Loading capacity        | 25 kg                                    |
| Laptop stand            | Up to 17"; Width: 424 mm; Length: 300 mm |

### Construction

|               |                 |
|---------------|-----------------|
| Main material | Aluminium, PMMA |
|---------------|-----------------|

### Color

|              |  |
|--------------|--|
| RAL 9016     |  |
| Custom color |  |

### Certificates

|                          |                             |
|--------------------------|-----------------------------|
| Certificates conforms to | IEC60601-1:2012 (Edition 3) |
| CE                       | Rohs compliant              |