



















# Budgie Range

Budgie Equipment Cart, Weight up to 13 Kgs.

MEDICAL  
CART COMPANY

Part Number	Image	Model Name	Comments	Datasheet	Packaging Dimension
BUG-0		Budgie Cart Basic	C2 Fix Height Column, 2 Locking Wheels 100mm, B1 Base		Box 1.: 620x450x180mm weight 8kgs Box 2.: 100x80x1000mm weight 3kgs
BUG-1		Budgie Cart with Vesa Mount	C2 Fix Height Column, 2 Locking Wheels 100mm, B1 Base, Display Mount		Box 1.: 620x450x180mm weight 8kgs Box 2.: 100x80x1400mm weight 4kgs
BUG-2		Budgie Cart with Tablet Mount	C2 Fix Height Column, 2 Locking Wheels 100mm, B1 Base, Tablet Mount		Box 1.: 620x450x180mm weight 8kgs Box 2.: 100x80x1400mm weight 4kgs
BUG-6		Budgie Cart with Quick Release Head	C2 Fix Height Column, 2 Locking Wheels 100mm, B1 Base, Quick Release Mount		Box 1.: 620x450x180mm weight 8kgs Box 2.: 100x80x1000mm weight 3kgs
BUG-7		Budgie Cart with Drip Pole	C2 Fix Height Column, 2 Locking Wheels 100mm, B1 Base, IV1 Drip Pole		Box 1.: 620x450x180mm weight 8kgs Box 2.: 100x80x1400mm weight 3kgs
BUG-2-US0		Budgie Ultrasound Cart	C2 Fix Height Column, 2 Locking Wheels 100mm, B1 Base, Tablet Mount, Gel Holder, Scanner Support		Box 1.: 620x450x180mm weight 9kgs Box 2.: 100x80x1400mm weight 3kgs
BUG-0-WS13		Small EKG, Fetal Monitor Budgie Cart	C2 Fix Height Column, 2 Locking Wheels 100mm, B1 Base, Blue Work Surface		Box 1.: 620x450x180mm weight 9kgs Box 2.: 100x80x1000mm weight 3kgs
BUG-0-WS12		Ultrasound Budgie Cart	C2 Fix Height Column, 2 Locking Wheels 100mm, B1 Base, Blue Work Surface (Rubber Scanner holders must be ordered separately)		Box 1.: 620x450x180mm weight 9kgs Box 2.: 100x80x1000mm weight 3kgs