

# Chameleon 4

Lightweight cart with sit-to-stand adjustment



## Ergonomic

Simple lightweight design for easy sit to stand adjustment



## Smooth and easy to clean surfaces

Quick and easy cleaning



## Durable

Lifetime Warranty\*




## Quick setup

Fast and simple to assemble



## Configuration options



**Display Mount** upto 24" Display



**Tablet Mount** 9 to 11" Tablets




**Laptop Mount** security mount 12-19"



□ Sit to Stand height adjustment

□ Sliding mount to adjust the height of your accessories



## Mounting options



**Display Mount**  
DM1



**Display Mount**  
DM3



**Display Mount 5**  
DM5



**Tablet Mount**  
TM1



**Tablet Mount**  
TM2



**Tablet Mount**  
CM1




**Tablet Mount**  
CM2



**Laptop Mount**  
LM1


## Keyboard Tray and accessories options



**Keyboard Tray**  
KB1



**Keyboard Tray**  
KB2



**2nd Work Surface**  
2WS4

## Holder options



**Front Din Rail**  
DR2, DR3, DR4, DR5, DR6



**Scanner holder**  
SH1,SH2



**Gel holder**  
GH1, GH2, GH3




**Fisher and Paylel**  
AM1



**Cable and Tube Hook**  
CH1

## Storage options




**Basket S**  
ST1



**Baskets L**  
ST2



**Storage Bin**  
ST3



**Cable Holder**  
ST4



**Handle and Storage**  
ST5



**CPU Mount Small**  
CPU1




**CPU Mount Medium**  
CPU2



**CPU Mount Big**  
CPU3



**Shelves**  
S1, S2, S3, S4  
S5, S6, S7



**PSU holder**  
PSH1, PSH2, PSH3,  
PSH4, PSH5, PSH6,  
PSH7

# Specifications



## Mechanical data

Height adjustment rate	Sit 775mm, stand 1135mm 360mm height adjustment
Worksurface	485x540mm
Wheelbase dimension	616mm (diameter)
Wheel construction	Polyamide with Polyurethane tread
Wheels/locking	4 lockable, 100mm diameter
Weight of a cart	18 kg
Loading capacity	25 kg

## Construction

Main material	Aluminum
Worksurface	Cast acrylic

## Color

RAL 9016	
Custom color	

## Certificates

Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant

## Technical characteristics

- Fully recycable
- Lifetime warranty\*
- Optional keyboard and mouse tray
- Large range of accessories
- Made in Hungary

