

Chameleon 4

Light weight cart with sit-to-stand adjustment



Ergonomic

Simple lightweight design for easy sit to stand adjustment



Smooth and easy to clean surfaces

To aid quicker and easier cleaning




Durable

Lifetime Warranty




Quick setup

Supplied fully assembled with minimal packaging



Configuration options






Display Mount upto 24" Display




Tablet Mount 9 to 11" Tablets



Laptop Mount security mount 12-19"





- ▣ Height adjustable monitor holder
- ▣ Modul Mount
- ▣ Handle for moving
- ▣ Sliding mount to adjust the height of your accessories

Laptop Security (12 to 19")

2nd Work surface - allow for full use of the work surface and easy access to the laptop

Keyboard and Mouse tray with full ergonomic adjustment

Range of storage solutions available (Baskets, bins)



Mounting options



Display Mount
DM1



Laptop Mount
LM1



Tablet Mount
TM1



Camera Mount
CM1



Module Mount
MM1

Keyboard options



2nd Worksurface
2WS4



Keyboard Tray
KB1

Din Rail options



Side Din Rail Kit
DR1



Scanner Holder
SH1



Scanner Holder
SH2




Mouse Holder
MH1




Gel Holder
GH1

Storage options



Basket S
ST1



Basket L
ST2



General Storage
ST4

Cabel, Tube and PSU Holder options



PSU Holder
PSH1



Hook
CH1

Wheels options



Antistatic options
100mm wheel



Mechanical data

Height adujstment rate	Sit 770mm, stand 1120mm 350mm height adjustment
Worksurface	500x660mm
Wheelbase dimension	616mm (diameter)
Wheel construction	Polyamide with Polyurethane tread
Wheels/locking	4 locking, 100mm diameter
Weight of a cart	18 kg
Loading capacity	25 kg

Construction

Main assembly	Aluminum
Worksurface	Cast acrylic

Color

RAL 9016	
Cutom color	

Mounting options

Display mount VESA 75/100
Laptop (secure) 12-19"
Tablet 9-11"
Camera Mount
Module Mount

Certificates

Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant

Technical characteristics

- Fully recycable
- Lifetime warranty
- Optional keyboard and mouse tray
- AC power pack option
- Large range of accessories
- Made in Hungary

