

Chameleon 5

Lightweight cart with sit-to-stand adjustment



Ergonomic

Simple lightweight design for easy sit to stand adjustment



Smooth and easy to clean surfaces

Quicker and easier cleaning



Multi-purpose

Easy integration in a variety of applications
(choice of laptop or VESA mount)



Durable

Lifetime Warranty*




Quick setup

Fast and simple to assemble



Configuration options



Display mount upto 24" Display




Tablet mount 9 to 11" Tablets




Laptop mount security mount 12-19"



Integrated Cable Holder



Dual purpose scanner Holder / Handle



Sit to Stand height adjustment



Height adjustable monitor Holder with built in tilt, swivel and rotate function

Sliding mount to adjust the height of your accessories




Optional Modul Mount

Optional storage for boxes, cleaning materials cables

Optional Din Rail to allow fitment of accessories

Range of storage solutions available (Baskets, bins, drawers)





- Height adjustable monitor Holder
- Optional Modul Mount
- Handle for moving
- Range of storage solutions available (Basket, bins, drawers)



Laptop Security (12 to 19")

2nd Work surface - allow for full use of the work surface and easy access to the laptop

Keyboard and Mouse tray with full ergonomic adjustment

Sliding mount of adjust the height of your accessories

Mounting options



Display Mount
DM1



Display Mount
DM3



Display Mount 5
DM5



Tablet Mount
TM1



Tablet Mount
TM2



Tablet Mount
CM1




Tablet Mount
CM2



Laptop Mount
LM1

Keyboard Tray and accessories options



Keyboard Tray
KB1



Keyboard Tray
KB2



2nd Work Surface
2WS4

Holder options



Front Din Rail
DR2, DR3, DR4, DR5, DR6



Scanner holder
SH1,SH2




Gel holder
GH1, GH2, GH3



Cable and Tube Hook
CH1

Storage options




Basket S
ST1



Baskets L
ST2



Storage Bin
ST3



Cable Holder
ST4



Handle and Storage
ST5



CPU Mount Small
CPU1




CPU Mount Medium
CPU2



CPU Mount Big
CPU3



Shelfes
S1, S2, S3, S4
S5, S6, S7



PSU holder
PSH1, PSH2, PSH3,
PSH4, PSH5, PSH6,
PSH7



Mechanical data

Height adjustment rate	Sit 780mm, stand 1140mm 360mm height adjustment
Worksurface	560x620mm
Wheelbase dimension	697mm (diameter)
Wheel construction	Polyamide with Polyurethane tread
Wheels/locking	5 lockable, 100mm diameter
Weight of a cart	21 kg
Loading capacity	25 kg

Construction

Main assembly	Aluminum
Worksurface	Cast acrylic

Color

RAL 9016	
Custom color	

Certificates

Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant

Technical characteristics

- Fully recycable
- Lifetime Warranty*
- Optional keyboard and mouse tray
- Large range of accessories
- Made in Hungary

