

Giraffe

Light weight cart



Ergonomic

Comfortable lightweight design



Smooth and easy to clean surfaces

To aid quicker and easier cleaning



Multi purpose

Easy integration in a variety of applications
(choice of tablet or VESA mount)



Durable

Lifetime Warranty



Quick setup

Supplied fully assembled with minimal packaging




Configuration options



Tablet Mount
9 to 11" Tablets



Top Plate



Quick Release Mount



Display Mount
upto 24" Display



Height adjustable monitor holder with built in tilt, swivel and rotate function

Handle for moving


Shelf (storage), a flat horizontal surface used

Range of storage solutions available

Height adjustable column

Basket

Custom color



Mounting options



Tablet Mount
TM1



Display Mount
DM1



Quick Release Mount
QR1



Top Plate with Din Rail
1911-9201

Keyboard options



Keyboard Tray
KB2

Din Rail Kit options



Front Din Rail
DR2, DR3, DR4, DR5



Scanner Holder
SH1



Mouse Holder
MH1



Gel Holder
GH1, GH2



Drip Pole
IV1



Din Rail Mount
DRC1, DRC2

Handle options



Handle
H1

Storage options



Basket S
ST1



Basket L
ST2



Storage Bin
ST3




Small Shelf
S7



Plastic Storage Bin
ST5

Cabel, Tube and PSU options



PSU Holder
PSH1



Hook
CH1

Wheels and Base options



Antistatic options
100mm wheels



Antistatic options
125mm wheels



Base options
B1



Base options
B3



Mechanical data

Height	1525mm
Wheelbase dimension	616mm (diameter)
Wheel construction	Polyamide with Polyurethane tread
Wheels/locking	4 locking, 100mm diameter
Loading capacity	20 kg
Angle of tilt	10°

Construction

Main assembly	Aluminium
---------------	-----------

Color

RAL 9016	
Custom color	

Mounting options

Display mount VESA 75/100

Tablet 9-11"

Top Plate

Quick Release Mount

Certificates

Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant

Technical characteristics

- Fully recyclable
- Lifetime Warranty
- AC power pack option
- Large range of accessories
- Made in Hungary

