

# Rhino

Equipment weight cart with sit-to-stand adjustment



## Ergonomic

Simple lightweight design for easy sit to stand adjustment



## Smooth and easy to clean surfaces

To aid quicker and easier cleaning



## Multi purpose

Easy integration in a variety of applications  
(choice of laptop or VESA mount)




## Durable

Lifetime Warranty




## Quick setup

Supplied fully assembled with minimal packaging



## Configuration options





**Display mount** upto 24" Display



**Tablet mount** 9 to 11" Tablets



**Laptop mount** security mount 12-19"



## Mounting options



Display Mount  
DM1



Laptop Mount  
LM1



Tablet Mount  
TM1



Camera Mount  
CM1

## Keyboard options



2nd Worksurface  
2WS4, 2WS5



Keyboard Tray  
KB1

## Din Rail options



Side Din Rail Kit  
DR1



Scanner holder  
SH1



Scanner holder  
SH2



Mouse holder  
MH1



Gel holder  
GH1, GH2

## Drawer options



Drawer



Medical Drawer



Pharmacy

## Storage options



Basket S  
ST1



Basket L  
ST2




Cabel holder  
ST4



Shelf  
S2, S3

## Cabel, Tube and PSU holder options



PSU holder  
PSH1



Hook  
CH1

## Wheels options



Antistatic options  
125mm wheel

## Wheels options



Achrily



Aluminium



## Mechanical data

Height adujstment rate	Sit 795mm, stand 1033mm 300x50mm height adjustment
Worksurface	557x638mm
Wheelbase dimension	749mm (diameter)
Wheel construction	Polyamide with Polyurethane tread
Wheels/locking	4 locking, 125mm diameter
Weight of a cart	21 kg
Loading capacity	60 kg

## Construction

Main assembly	Aluminum
Worksurface	Cast acrylic

## Color

RAL 9016	
Cutom color	

## Mounting options

Display mount VESA 75/100
Laptop (secure) 12-19"
Tablet 9-11"
Camera Mount

## Certificates

Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant

## Technical characteristics

- Fully recycable
- Lifetime Warranty
- Optional keyboard and mouse tray
- AC power pack option
- Large range of accessories
- Made in Hungary

