

Chamelen 4 Cart Ultrasound

MEDICAL
CART COMPANY

CH4-0-WS12



Ergonomic

Simple, light weight, easy to use



Smooth and easy to clean surfaces

Quick and easy cleaning



Multi-purpose

Customise with a choice of accessories




Durable

Lifetime Warranty*



Quick setup


Supplied fully assembled with minimal packaging



Gel Holder




9-11" Tablet Mount




Scanner Holder





Lock with key



Optional Storage



Cable Management



Height Adjustable



Mechanical data

Height	Min: 760mm Max: 1110mm
Wheelbase dimension	616mm (diameter)
Wheel construction	Polyamide with Polyurethane tread
Wheels/locking	4 lockable, 100mm diameter
Loading capacity	20 kg

Construction

Main assembly	Aluminium
---------------	-----------

Color

RAL 9016	
Custom color	

Accessories options

ST3 Storage Bin

ST2 Basket

CH1 Cable Hook

Certificates

Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant

Technical characteristics

- Fully recyclable
- Lifetime Warranty*
- Made in Hungary

