

Chameleon 4 Cart DM5 Display Mount

CH4-DM5

MEDICAL
CART COMPANY



Ergonomic

Simple, light weight, easy to use



Smooth and easy to clean surfaces

Quick and easy cleaning



Multi-purpose

Customise with a choice of accessories




Durable


Lifetime Warranty*



Quick setup

Supplied fully assembled with minimal packaging






Adjustable Monitor height

- Max: 414mm
- Min: 164mm



Adjustable Column height

- Max Work Surface height: 1135mm
- Min Work Surface height: 775mm



Specifications

Mechanical data

Height	Min: 1175mm Max: 1540mm
Wheelbase dimension	616mm (diameter)
Wheel construction	Polyamide with Polyurethane tread
Wheels/locking	4 lockable, 100mm diameter
Loading capacity	20 kg

Construction

Main assembly	Aluminium
---------------	-----------

Color

RAL 9016	
Custom color	

Certificates

Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant

Technical characteristics

- Fully recyclable
- Lifetime Warranty*
- Made in Hungary

