

Two Person Work Station

CHA5-0-TS

MEDICAL
CART COMPANY



Ergonomic



Quick and easy to clean



Multi-purpose



Lifetime Warranty



Quick setup



Mechanical data

Height	1135mm
Wheelbase dimension	616mm (diameter)
Wheel construction	Polyamide with Polyurethane tread
Wheels/locking	2 lockable, 697mm diameter
Loading capacity	20 kg

Construction

Main material	Aluminum
---------------	----------

Color

RAL 9016	
Custom color	

Certificates

Certificates conforms to	IEC60601-1:2012 (Edition 3)
CE	Rohs compliant

Technical characteristics

- Fully recyclable
- Lifetime Warranty
- Made in Hungary

