VENTILATOR CART SOLUTIONS



Medical Cart
COMPANY

GENERAL FEATURES

EASY TO CLEAN SURFACES
Clean surfaces without dirt traps.
One wipe cleaning.

QUICK SETUP
Available assembled or
flat-packed, easy to assemble.

7 YEARS WARRANTY
Seven-year warranty for
the care-free treatment.

ERGONOMIC
Designed for easy
and comfortable use.

SMALL FOOTPRINT
Developed to fit narrow doors,
ease of movement.


LOCKABLE WHEELS
Simple locking mechanism
operated by foot.





CUSTOMIZED WORKTOPS

We customize the worktop design
to the ventilator product used on it.

If our customers provide the mounting data
of the device, our engineering team will fit
the top plate perfectly.







Visit the Tiger webshop and select
the color you desire:
https://www.tiger-coatings.com


Datasheet
Scan the QR

CART EXAMPLES

VENTILATOR 1 custom cart

Custom design sheet metal worktop with DIN rails +

- Handle
- Accessory mount
- Cable hook



VENTILATOR 1 custom cart

Custom design sheet metal worktop with DIN rails +

- Handle with hook
- Accessory mount (2x)
- Drip pole
- Basket



VENTILATOR 2 custom cart

Custom design worktop with wrap-around DIN rail frame +

- Snake 1 arm
- Accessory mount
- Basket
- Cable hook



Images shown here are for illustrative purpose only.

CART EXAMPLES

VENTILATOR 3 custom cart



Custom design worktop
with wrap-around
DIN rail frame +

- Snake 1 arm
- Accessory mount
- Gas bottle mount

VENTILATOR 4 custom cart



Custom design worktop
with wrap-around
DIN rail frame +

- Snake 1 arm
- Accessory mount
- Cable hook

CORE PRODUCTS

VEN1-0

BASIC LIGHTWEIGHT CART

- Lightweight base
- 100mm diameter wheels
- Lightweight column
- Top mount plate according to customer's need (see page 5)
- Accessory options (see page 7)

Datasheet
Scan the QR



VEN1-7

LIGHTWEIGHT GAS BLENDER CART

- Lightweight base
- 100mm diameter wheels
- Lightweight column
- Drip pole inside the column
- Accessory options (see page 7)

Datasheet
Scan the QR




VEN2-0

BASIC MIDWEIGHT CART

- Midweight base
- 100mm diameter wheels
- Strong column
- Top mount plate according to customer's need (see page x)
- Accessory options (see page x)

Datasheet
Scan the QR



CORE PRODUCTS

VEN3-0

BASIC MIDWEIGHT CART WITH GAS BOTTLE HOLDER

- Strong base with 2x gas bottle holder
- 125mm diameter wheels
- Strong column
- Top mount plate according to customer's need (see page 5)
- Accessory options (see page 7)

Datasheet
Scan the QR



VEN4-0

BASIC COMPRESSOR CART WITH GAS BOTTLE HOLDER

- Strong base with Compressor and 2x gas bottle holder
- 100mm diameter wheels
- Strong column
- Top mount plate according to customer's need (see page 5)
- Accessory options (see page 7)


Datasheet
Scan the QR




QUICK RELEASE MOUNT

We offer a series of Quick release mounts for our carts.
The plate that fits the mount can be customized.
Any device's mounting hole pattern can be added to the plate.

Basic mounts
with spring locking pin





Premium mount
with spring loaded release button




Datasheet
Scan the QR

MOUNTING ACCESSORIES

The columns are designed for fitting accessories on the front and back. Channels run all along the column. Our accessories can be mounted on the column through sliding nuts with M6 screws.






The Top plates have DIN rails on the sides. We offer a series of DIN rail clamps.

Datasheet
Scan the QR

CIRCUIT ARM

- Free movement with 5 joints
- Lightweight structure
- Several accessory mounting options
- Cable/tube leading clip-in options



MOUNTED
ON DIN RAIL

TUBE HOLDER
CLIP-IN

FLEXIBLE END
MOUNT WITH BALL
JOINT TUBE HOLDER
- MAX. LOAD: 0.9kg

Datasheet
Scan the QR

ACCESSORIES


CIRCUIT ARM



SNA1-WS2-L2B-TM2

Datasheet
Scan the QR

WORKTOP



CUSTOM DESIGN

Datasheet
Scan the QR

STORAGE



ST1/2/3/5/6/7/8/9

Datasheet
Scan the QR

DIN RAIL



DR2, DR3

Datasheet
Scan the QR

ACCESSORY MOUNT



AM1-C5

Datasheet
Scan the QR

GAS BOTTLE HOLDER



GBH 1/2

Datasheet
Scan the QR

DIN RAIL CLAMP



DRC 1/2/3/4

Datasheet
Scan the QR

CABLE HOOK



CH1, CH2

Datasheet
Scan the QR

DRAWERS



D1/2/3/5

Datasheet
Scan the QR

QUICK RELEASE MOUNT



QRM

Datasheet
Scan the QR

Medical Cart COMPANY



Medical Cart Company Kft.
2051 Hungary
Biatorbágy, Budai út 16.
Phone: +36 23 920 790
Email: info@medicalcartcompany.hu
www.medicalcartcompany.hu

Revision00 - 2026.10.12.